


iMEDICOM
Intended for Human



EZ Sinuscare

Balloon sinuplasty system



iMEDICOM
Intended for Human

www.imedicom.co.kr


What is a Balloon Sinuplasty?

A balloon sinuplasty is a minimally invasive endoscopic procedure used to treat chronic sinusitis.

Using a small, flexible balloon catheter to open up blocked sinus passageways to relieve pressure and pain in patients suffering from chronic sinusitis symptoms.

Product Overview

Sinuscore EZ



○ Balloon Catheter



Balloon length	Balloon diameter	Balloon pressure
10, 20mm	5~9mm	12ATM (176psi)

○ Guide

Type	Degree	
Sphenoid	0°	
Frontal	70°	
Maxillary	110°	

Indications

- Chronic rhinosinusitis
- Subacute rhinosinusitis
- When medication is not effective
- Revision surgery
- Inflammatory and infectious sinus diseases

Advantages

- Absence of surgical incisions
- Lack of serious complications
- Safe and effective
- Not necessary to remove any bone and tissue
- Short procedure (15mins to 30mins)
- Return to normal activities within 24 hours

Procedure

Prepare an Inflation device and fill it up with sterile saline.


Connect the end of the Inflation device to the Y-port of the Balloon catheter.

Choose an appropriate Guide for each targeted sinus (0° : Sphenoid, 70° : Frontal, 110° : Maxillary).

Attach the Guide to the outer tube of the Balloon catheter.

Push the Position handle back for the balloon to be located in the Guide.

01




Insert the tube along the blocked sinus by gentle movements.
(Do not use excessive force.)

Advance the Balloon catheter over the Guide and visually confirm that the Balloon has entirely exited the Guide through endoscope.

Inflate the balloon up to 12 atmospheres with a recommended total hold time of 20 seconds.

* Do not exceed 12ATM (Max. pressure).

02



Deflate the Balloon catheter with Inflation device.

Push the Position handle back for the balloon to be located in the Guide.

Remove the Balloon catheter from the patient.

* Ensure if the balloon is completely deflated before moving the balloon catheter. Advancing or retracting the sinus balloon while it is inflated may cause serious damage to anatomy or to the product.

Ordering Information

Sinuscore EZ instrument

Product	Order code	Description		
		Balloon length	Balloon diameter	Guide
Sinuscore EZ	SCLZ105	10mm	5mm	0° 70° 110°
	SCLZ106		6mm	
	SCLZ107		7mm	
	SCLZ108		8mm	
	SCLZ109		9mm	
	SCLZ205	20mm	5mm	
	SCLZ206		6mm	
	SCLZ207		7mm	
	SCLZ208		8mm	
	SCLZ209		9mm	

Accessory (to be purchased separately)

Product	Order code	Description
Inflation device	IBE-B00	0 ATM to 24.6 ATM



IMEDICOM CO., LTD.

#612, LS-ro 172, Gunpo-si, Gyeonggi-do, 15807 Republic of Korea
 Tel: +82 31 479 1156 Fax: +82 31 479 1157
 www.imedicom.co.kr

Distributed by