









Ref. No.	Description	
K6T. F90. 00	PTS Tray	
K6T. F91. 00	Bone Screw Tray	
K6T. F30	Impactor	
CK6T. F50. 11 CK6T. F50. 18	Drill Bit Φ1.1 Drill Bit Φ1.8	
K6T. F60. 15 K6T. F60. 27 K6T. F60. 01	Screwdriver Φ1.5 Screwdriver Φ2.7 Screwdriver Handle	
K6T. F70. 15 K6T. F70. 27	Depth Gauge Φ1.5 Depth Gauge Φ2.7	
K6T. F80	Drill Guide	
K6T. F10. 04	Trial #4	
K6T. F10. 05	Trial #5	
K6T. F10. 06	Trial #6	
K6T. F10. 07	Trial #7	
K6T. F10. 08	Trial #8	
K6T. F10. 09	Trial #9	
K6T. F10. 10	Trial #10	
K6T. F10. 11	Trial #11	
K6T. F10. 12	Trial #12	
K6T. F93. 01 K6T. F93. 02	Templates 100%	
K6T. F93. 03 K6T. F93. 04	Templates 105%	
K6T. F93. 05 K6T. F93. 06	Templates 110%	


Aniverse™ Patella Tracking System

Many pet dogs suffer from habitual patellar luxation and dislocation. Aniverse™ Patella Tracking System is veterinary medical device designed to help you enjoy healthy walks with your canine companions.



DOC NO. : VTMSC-0101-M04

Aniverse™ Patella Tracking System



01 Trochlear Cap PEEK
(Polyetheretherketon : ASTM F2026)
> Biocompatible product

02 Screw
Ti-6AL-4V ELI (ASTM F136)
> Self tapping screw


03 Baseplate
Ti-6AL-4V ELI (ASTM F136)
> Surface treatment improve optimal osseointegration

Patella Tracking System



Trochlear Cap_Ref. No.	For use with	Baseplate_Ref. No.
K6I. F10. 04	Φ1.5 Bone Screws	K6I. F20. 04
K6I. F10. 05		K6I. F20. 05
K6I. F10. 06		K6I. F20. 06
K6I. F10. 07		K6I. F20. 07
K6I. F10. 08	Φ2.7 Bone Screws	K6I. F20. 08
K6I. F10. 09		K6I. F20. 09
K6I. F10. 10		K6I. F20. 10
K6I. F10. 11		K6I. F20. 11
K6I. F10. 12		K6I. F20. 12

Bone Screws



Φ1.5 Bone Screws for use with PTS No. 4~7		Φ2.7 Bone Screws for use with PTS No. 8~12	
Ref. No.	Description	Ref. No.	Description
K6I. F31. 06	6mm	K6I. F32. 10	10mm
K6I. F31. 07	7mm	K6I. F32. 12	12mm
K6I. F31. 08	8mm	K6I. F32. 14	14mm
K6I. F31. 10	10mm	K6I. F32. 16	16mm
K6I. F31. 12	12mm	K6I. F32. 18	18mm
K6I. F31. 14	14mm	K6I. F32. 20	20mm
K6I. F31. 16	16mm	K6I. F32. 22	22mm
K6I. F31. 18	18mm	K6I. F32. 24	24mm
K6I. F31. 20	20mm	K6I. F32. 26	26mm
K6I. F31. 22	22mm	K6I. F32. 28	28mm
K6I. F31. 24	24mm		